



## Nickel acid baths cuvette test 5-120 g/L Ni, 25 tests

Product #: LCK237

AED Price: Contact Hach



## LCK cuvette tests: Minimum amounts of reagents, maximum safety.

Pre-dosed reagents for maximum safety. Easy to use and photometric evaluation, packaged in a fully equipped box. Environmentally friendly due to less chemicals.

## **Specifications**

Description: Acidic Nickel plating baths, Nickel sulphate and nickel sulphamate baths

Instrument: DR3900, DR6000, DR1900, DR2800, DR3800, DR5000

Measuring range: 5 - 120 g/L Ni

Method: Intrinsic Baths Colour

Number of tests: 25

Parameter: Nickel Platform: LCK

Storage conditions: 15 - 25 °C

## **Required Accessories**

• Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)