



Zirconium cuvette test, 6-60 mg/L Zr, 25 tests

Product #:

LCK364

AED Price:

Contact Hach

Hazardous



LCK cuvette tests: Minimum amounts of reagents, maximum safety.

The first specific, easy to use rapid test which helps to evaluate the Zirconium concentration during the passivation process.

Pre-dosed reagents for maximum safety. Easy to use and photometric evaluation, packaged in a fully equipped box. Environmentally friendly due to less chemicals.

Helps adjusting the right amount of Zirconium within the process

Helps minimising chemistry consumption and costs in Zirconium based coating processes

Helps ensuring a consistent high product quality

Specifications

| | |
|---------------------|--|
| EPA compliant: | N/A |
| Instrument: | DR3900, DR6000, DR1900, DR2800, DR3800, DR5000 |
| Measuring range: | 6 - 60 mg/L |
| Method: | SurTec/Hach Method |
| Number of tests: | 25 |
| Parameter: | Zirconium |
| Platform: | LCK |
| Storage conditions: | 2 - 8 °C (keep refrigerated) |

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)