



Free Chlorine cuvette test, 0.05-2.0 mg/L Cl_2 , 24 tests

Product #:

LCK410

AED Price:

Contact Hach

Hazardous



LCK cuvette tests: Minimum amounts of reagents, maximum safety.

Cuvette test for the determination of free chlorine, 0.05-2.0 mg/L Cl_2 and chlorine dioxide, 0.09-3.8 mg/L ClO_2 .

Pre-dosed reagents for maximum safety. Easy to use and photometric evaluation, packaged in a fully equipped box. Environmentally friendly due to less chemicals.

Precise and reliable measurement values

Maximum safety for users, thanks to the closed system and low amounts of reagents

Convenient and error-free dosing of the reagents

Barcode label for automatic recognition in the photometer

Differentiation between tests and measuring ranges by means of colour-coding

Specifications

According to standard:	ISO 7393-1-2-1985, DIN 38408 G4-2
Instrument:	DR3900, DR6000, DR1900, DR2800, DR3800, DR5000
Measuring range:	0.05 - 2.0 mg/L Chlorine free/ ClO_2
Method:	DPD
Number of tests:	24 (test includes zero solution)
Parameter:	Chlorine, free

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)