



Silica acid (high range) reagent set 0.8-100 mg/L SiO₂

Product #:

LCS028

AED Price:

Contact Hach

Hazardous



Or 0.4-50 mg/L Si

Specifications

Description:	reagent set for Silica acid/Silicon high range
Measuring range:	0.8 - 100 mg/L SiO ₂
Measuring range (2):	0.4 - 50 mg/L Si
Method:	Molybdenum blue
Number of tests:	50
Parameter:	Silica acid
Sample size:	5 mL
Special Features:	HR, 50 tests
Special remarks:	round cuvette LCW906 needed
Storage conditions:	15 °C - 25 °C

Required Accessories

- Round cuvettes, glass, 13 mm, 25 pcs (Item LCW906)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)