



Acid capacity - KS4.3, cuvette test, 0.5-8.0 mmol/L, 25 tests

Product #: LCK362
AED Price: Contact Hach

LCK cuvette tests: Minimum amounts of reagents, maximum safety.

Ready for use, ensures maximum safety, easy to handle, photometric evaluation system and environmentally friendly

Precise and reliable measurement values

Maximum safety for users, thanks to the closed system and low amounts of reagents

Convenient and error-free dosing of the reagents

Barcode label for automatic recognition in the photometer

Differentiation between tests and measuring ranges by means of colour-coding

Specifications

Description:	KS _{4.3} /Calcium Carbonate CaCO ₃
EPA compliant:	N/A
Instrument:	DR3900, DR6000, DR1900, DR2800, DR3800, DR5000
Measuring range:	0.5 - 8.0 mmol/L
Measuring range (2):	25 - 400 mg/L CaCO ₃
Method:	Hach Method
Number of tests:	25
Parameter:	Acid capacity
Platform:	LCK
Storage conditions:	15 - 25 °C

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)