



COD cuvette test 1,000-10,000 mg/L O₂, 25 tests

Product #:

LCK014

AED Price:

Contact Hach

Hazardous



LCK cuvette tests: Minimum amounts of reagents, maximum safety.

Pre-dosed reagents for maximum safety. Easy to use and photometric evaluation, packaged in a fully equipped box. Officially approved for consent limits. Environmentally friendly due to less chemicals.

Precise and reliable measurement values

Maximum safety for users, thanks to the closed system and low amounts of reagents

Convenient and error-free dosing of the reagents

Barcode label for automatic recognition in the photometer and Truecal calibration

Differentiation between tests and measuring ranges by means of colour-coding

Specifications

According to standard:	ISO 6060-1989, DIN 38409-H41-H44
Description:	Chemical Oxygen Demand
Digestion required:	Yes
Instrument:	DR3900, DR6000, DR1900, DR2800, DR3800, DR5000
Measuring range:	1000 - 10000 mg/L O ₂
Method:	Dichromate
Number of tests:	25
Parameter:	COD
Platform:	LCK
Storage conditions:	15 - 25 °C (protect from light)
Truecal:	Yes

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)