



## BOD5 cuvette test 4 -1650 mg/L O<sub>2</sub>, 39 tests

Product #:

LCK555

AED Price:

Contact Hach

Hazardous



### LCK cuvette tests: Minimum amounts of reagents, maximum safety.

Pre-dosed reagents for maximum safety. Easy to use and photometric evaluation, packaged in a fully equipped box. Environmentally friendly due to less chemicals.

#### Precise and reliable measurement values

Maximum safety for users, thanks to the closed system and low amounts of reagents

Convenient and error-free dosing of the reagents

Barcode label for automatic recognition in the photometer

Differentiation between tests and measuring ranges by means of colour-coding

---

### Specifications

According to standard:	EN 1899-1
Description:	Biological Oxygen Demand, 5 days
EPA compliant:	N/A
Instrument:	DR3900, DR6000, DR1900, DR2800, DR3800, DR5000
Measuring range:	4 - 1650 mg/L O <sub>2</sub>
Method:	Dilution Method
Number of tests:	39
Parameter:	BOD <sub>5</sub>
Platform:	LCK
Storage conditions:	2 - 8 °C (keep refrigerated)

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)