



Starch cuvette test 2-150 mg/L, 25 tests

Product #:

LCK357

Hazardous

AED Price:

Contact Hach



Pre-dosed reagents for maximum safety. Easy to use and photometric evaluation, packaged in a fully equipped box. Environmentally friendly due to less chemicals.

Specifications

Description:	cuvette test for soluble Starch
EPA compliant:	N/A
Instrument:	DR3900, DR6000, DR1900, DR2800, DR3800, DR5000
Measuring range:	2 - 150 mg/L Starch
Method:	Hach Method
Number of tests:	25
Parameter:	Starch
Platform:	LCK
Storage conditions:	15 - 25 °C

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)