



Carbonate/carbon dioxide cuvette test 55-550 mg/L CO₂, 25 tests

Product #:

LCK388

AED Price:

Contact Hach

LCK cuvette tests: Minimum amounts of reagents, maximum safety.

Ready for use, for maximum safety, easy handling and photometric evaluation, environmentally friendly

Precise and reliable measurement values

Maximum safety for users, thanks to the closed system and low amounts of reagents

Convenient and error-free dosing of the reagents

Barcode label for automatic recognition in the photometer

Differentiation between tests and measuring ranges by means of colour-coding

Specifications

Digestion required:	yes
EPA compliant:	N/A
Instrument:	DR3900, DR6000, DR1900, DR2800, DR3800, DR5000
Measuring range:	55 - 550 mg/L CO ₂
Method:	pH Indicator
Number of tests:	25
Parameter:	Carbonate/carbon dioxide
Platform:	LCK
Storage conditions:	15 - 25 °C

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)