



Sulphide cuvette test 0.1-2.0 mg/L S²⁻, 25 tests

Product #:

LCK653

AED Price:

Contact Hach

Hazardous



LCK cuvette tests: Minimum amounts of reagents, maximum safety.

Pre-dosed reagents for maximum safety. Easy to use and photometric evaluation, packaged in a fully equipped box. Environmentally friendly due to less chemicals.

Precise and reliable measurement values

Maximum safety for users, thanks to the closed system and low amounts of reagents

Convenient and error-free dosing of the reagents

Barcode label for automatic recognition in the photometer

Differentiation between tests and measuring ranges by means of colour-coding

Specifications

According to standard:	ISO 10530-1991, DIN 38405-D26
Description:	dissolved Sulphide
Instrument:	DR3900, DR6000, DR1900, DR2800, DR3800, DR5000
Measuring range:	0.1 - 2.0 mg/L S ²⁻
Method:	Dimethyl-p-phenylenediamine
Number of tests:	25
Parameter:	Sulphide
Platform:	LCK
Storage conditions:	15 - 25 °C

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)