



## Sulphite cuvette test 0.1-5.0 mg/L SO<sub>3</sub>, 25 tests

**Product #:** LCK654  
**AED Price:** Contact Hach

### LCK cuvette tests: Minimum amounts of reagents, maximum safety.

Pre-dosed reagents for maximum safety. Easy to use and photometric evaluation, packaged in a fully equipped box. Environmentally friendly due to less chemicals.

#### Precise and reliable measurement values

#### Maximum safety for users, thanks to the closed system and low amounts of reagents

#### Convenient and error-free dosing of the reagents

#### Barcode label for automatic recognition in the photometer

#### Differentiation between tests and measuring ranges by means of colour-coding

---

### Specifications

|                     |  |
|---------------------|--|
| EPA compliant:      | N/A  |
| Instrument:         | DR3900, DR6000, DR1900, DR2800, DR3800, DR5000 |
| Measuring range:    | 0.1 - 5.0 mg/L SO <sub>3</sub>                 |
| Method:             | Hach Method                                    |
| Number of tests:    | 25   |
| Parameter:          | Sulphite                                       |
| Platform:           | LCK  |
| Storage conditions: | 15 - 25 °C                                     |

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)