



Addista Multi parameter standard solution

Product #: LCA700
AED Price: Contact Hach

System for Analytical Quality Assurance (AQA), consisting of one combined standard/spiking solution and two round robin test solutions.

Applicable for the cuvette tests

- LCK304 Ammonium, 0.015-2.0 mg/L $\text{NH}_4\text{-N}$
- LCK311 Chloride, 1-70 mg/L Cl
- LCK228 Potassium, 5-50 mg/L K
- LCK328 Potassium, 8-50 mg/L K
- LCK348 Phosphate (ortho), 0.5-5.0 mg/L $\text{PO}_4\text{-P}$
- LCK414 / LCK1414 COD, 5-60 mg/L O_2
- LCK238 Laton, 5-40 mg/L TN_b

Specifications

Description:	ADDISTA Multi parameter standard solution
Platform:	LCA
Storage conditions:	2 - 8 °C (keep refrigerated)