



Addista Multi parameter standard solution

Product #: LCA701
AED Price: Contact Hach

System for Analytical Quality Assurance (AQA), consisting of one combined standard/spiking solution and two round robin test solutions.

Applicable for the cuvette tests

- LCK306 Lead, 0.1-2.0 mg/L Pb
- LCK321 Iron, 0.2-6.0 mg/L Fe
- LCK329 Copper, 0.1-8.0 mg/L Cu
- LCK337 Nickel, 0.1-6.0 mg/L Ni
- LCK353 Sulphate, 150-900 mg/L SO₄
- LCK360 Zinc, 0.2-6.0 mg/L Zn

Specifications

Description:	ADDISTA Multi parameter standard solution
Platform:	LCA
Storage conditions:	2 - 8 °C (keep refrigerated)