



Addista Multi parameter standard solution

Product #: LCA702
AED Price: Contact Hach

System for Analytical Quality Assurance (AQA), consisting of one combined standard/spiking solution and two round robin test solutions.

Applicable for the cuvette tests

- LCK301 Aluminium, 0.02-0.5 mg/L Al
- LCK308 Cadmium, 0.02-0.3 mg/L Cd
- LCK313 Chromium (VI), 0.03-1.0 mg/L Cr
- LCK313 Chromium (total), 0.03-1.0 mg/L Cr
- LCS313 Chromium trace, 0.005-0.25 mg/L Cr
- LCK353 Sulphate, 150-900 mg/L SO₄

Specifications

Description:	ADDISTA Multi parameter standard solution
Platform:	LCA
Storage conditions:	2 - 8 °C (keep refrigerated)