



Addista Multi parameter standard solution

Product #: LCA703
AED Price: Contact Hach

System for Analytical Quality Assurance (AQA), consisting of one combined standard/spiking solution and two round robin test solutions.

Applicable for the cuvette tests

- LCK049 Orthophosphate, 1.6-30 mg/L $\text{PO}_4\text{-P}$
- LCK114 COD, 150-1000 mg/L O_2
- LCI400 COD, 0-1000 mg/L O_2
- LCK303 Ammonium, 2-47 mg/L $\text{NH}_4\text{-N}$
- LCK311 Chloride, 1-70 mg/L Cl
- LCK339 Nitrate, 0.23-13.5 mg/L $\text{NO}_3\text{-N}$
- LCK350 Phosphate (ortho), 2-20 mg/L $\text{PO}_4\text{-P}$
- LCK353 Sulphate, 150-900 mg/L SO_4
- LCK386 TOC, 30-300 mg/L C

Specifications

Description:	ADDISTA Multi parameter standard solution
Platform:	LCA
Storage conditions:	2 - 8 °C (keep refrigerated)