



Addista Multi parameter standard solution

Product #:	LCA705
AED Price:	Contact Hach

System for Analytical Quality Assurance (AQA), consisting of one combined standard/spiking solution and two round robin test solutions.

Applicable for the cuvette tests

- LCK014 COD, 1000-10000 mg/L O₂
- LCK302 Ammonium, 47-130 mg/L NH₄-N
- LCK311 Chloride, 1-70 mg/L Cl
- LCK387 TOC, 300-3000 mg/L C

Specifications

Description:	ADDISTA Multi parameter standard solution
Platform:	LCA
Storage conditions:	2 - 8 °C (keep refrigerated)