



Addista Multi parameter standard solution

Product #: LCA708
AED Price: Contact Hach

System for Analytical Quality Assurance (AQA), consisting of one combined standard/spiking solution and two round robin test solutions.

Applicable for the cuvette tests

- LCK338 LATON, 20-100 mg/L TN_b
- LCK514 COD, 100-2000 mg/L O_2
- LCK350 Phosphate (total), 2-20 mg/L $\text{PO}_4\text{-P}$

Specifications

Description:	ADDISTA Multi parameter standard solution
Platform:	LCA
Storage conditions:	2 - 8 °C (keep refrigerated)