



## Addista Multi parameter standard solution

**Product #:** LCA709  
**AED Price:** Contact Hach

System for Analytical Quality Assurance (AQA), consisting of one combined standard/spiking solution and two round robin test solutions.

Applicable for the cuvette tests

- LCK138 LATON, 1-16 mg/L  $\text{TN}_b$
- LCK614 COD, 50-300 mg/L  $\text{O}_2$
- LCK349 Phosphate (total), 0.05-1.5 mg/L  $\text{PO}_4\text{-P}$
- LCK342 Nitrite, 0.6-6.0 mg/L  $\text{NO}_2\text{-N}$

---

### Specifications

|                     |                                           |
|---------------------|-------------------------------------------|
| Description:        | ADDISTA Multi parameter standard solution |
| Platform:           | LCA                                       |
| Storage conditions: | 2 - 8 °C (keep refrigerated)              |