



## Ammonium cuvette test 0.015-2.0 mg/L $\text{NH}_4\text{-N}$ , 25 tests

Product #:

LCK304

AED Price:

Contact Hach

Hazardous



### LCK cuvette tests: Minimum amounts of reagents, maximum safety.

Pre-dosed reagents for maximum safety. Easy to use and photometric evaluation, packaged in a fully equipped box. Officially approved for consent limits. Environmentally friendly due to less chemicals.

#### Precise and reliable measurement values

Maximum safety for users, thanks to the closed system and low amounts of reagents

Convenient and error-free dosing of the reagents

Barcode label for automatic recognition in the photometer and Truecal calibration

Differentiation between tests and measuring ranges by means of colour-coding

---

### Specifications

According to standard:	ISO 7150-1, DIN 38406 E5-1, UNI 11669:2017, ISO 23695
Instrument:	DR3900, DR6000, DR1900, DR2800, DR3800, DR5000
Measuring range:	0.015 - 2.0 mg/L $\text{NH}_4\text{-N}$
Measuring range (2):	0.02 - 2.5 mg/L $\text{NH}_4$
Method:	Indophenol Blue
Number of tests:	25
Parameter:	Ammonium
Platform:	LCK
Storage conditions:	2 - 8 °C (keep refrigerated)
Truecal:	Yes

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)