



Laton Total Nitrogen cuvette test 20-100 mg/L TN, 25 tests

Product #:

LCK338

AED Price:

Contact Hach

Hazardous



LCK cuvette tests: Minimum amounts of reagents, maximum safety.

Pre-dosed reagents for maximum safety. Easy to use and photometric evaluation, packaged in a fully equipped box. Officially approved for consent limits. Environmentally friendly due to less chemicals.

Precise and reliable measurement values

Maximum safety for users, thanks to the closed system and low amounts of reagents

Convenient and error-free dosing of the reagents

Barcode label for automatic recognition in the photometer and Truecal calibration

Differentiation between tests and measuring ranges by means of colour-coding

Specifications

According to standard:	EN ISO 11905-1, ISO23697-1
Description:	Total Nitrogen bonded
Digestion required:	Yes
EPA compliant:	N/A
Instrument:	DR3900, DR6000, DR1900, DR2800, DR3800, DR5000
Measuring range:	20 - 100 mg/L TN _b
Method:	Koroleff Digestion (Peroxodisulphate), and Photometric Detection with 2,6-Dimethylphenol
Number of tests:	25
Parameter:	Nitrogen total (Laton)
Platform:	LCK
Standard method:	EN ISO 11905-1 Digestion with Peroxodisulphate
Storage conditions:	15 - 25 °C
Truecal:	Yes

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)