



## Fluoride cuvette test 0.1-2.5 mg/L F, 25 tests

**Product #:**

**LCK323**

**Hazardous**

**AED Price:**

Contact Hach



### LCK cuvette tests: Minimum amounts of reagents, maximum safety.

Ready for use, ensures maximum safety, easy to handle, photometric evaluation system and environmentally friendly

**Precise and reliable measurement values**

**Maximum safety for users, thanks to the closed system and low amounts of reagents**

**Convenient and error-free dosing of the reagents**

**Barcode label for automatic recognition in the photometer**

**Differentiation between tests and measuring ranges by means of colour-coding**

---

### Specifications

EPA compliant:	N/A
Instrument:	DR3900, DR6000, DR1900, DR2800, DR3800, DR5000
Measuring range:	0.1 - 2.5 mg/L F
Method:	SPADNS
Number of tests:	25
Parameter:	Fluoride
Platform:	LCK
Storage conditions:	15 - 25 °C

---

### Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)