



Lumistox Luminescent bacteria test

Product #:	LCK482
AED Price:	Contact Hach

Luminescent bacteria for 400 toxicity tests, conserved in 20 small tubes, 250 mL storage bottle reactivation solution

Specifications

Description:	Luminescent bacteria for 400 toxicity tests, preserved in 20 small tubes, 250 ml storage bottle reactivation solution
Method:	Toxicity
Parameter:	LUMISTOX luminescent bacteria test
Standard method:	EN ISO 11348 part 2
Storage conditions:	-18 °C