



Lumismini Luminescent bacteria test

Product #:	LCK487
AED Price:	Contact Hach

Luminescent bacteria for 90 toxicity tests, preserved in 10 small tubes, 50 ml stock bottle reactivation solution

Specifications

Description:	Luminescent bacteria for 90 toxicity tests, preserved in 10 small tubes, 50 ml stock bottle reactivation solution
Method:	Toxicity
Parameter:	LUMISMINI liminiscent bacteria test
Standard method:	EN ISO 11348 part 2
Storage conditions:	-18 °C