



## DREL Industrial water lab

Product #: LZV735

**OBSELETE ITEM**  
This item is no longer available.

Portable water lab including spectrophotometer DR1900, complete with batteries and cuvettes, instrument case, reagent & apparatus case, reagent set, and instrument set.

Today, we have to be able to carry out analysis anywhere. The portable laboratory contains everything needed to ensure that accurate measurements can be taken.

Even in the field, there is the need for efficient measurement technology, practical accessories and ready-to-use reagents. The integrated work surface enables analysis to be carried out safely and easily at any location.

**Complete water lab for field testing, incl. photometer, digital titrator and electrochemical methods**

**Versatile and rugged**

**Simple operation with ready-to-use reagents**

---

### Specifications

Height:	280 mm
Includes:	VIS spectrophotometer DR 2800 (LPV422.99.00001), instrument case (5838500), reagent case (5839400), reagent set (232330), reagent set (2923400)
Wavelength range:	340 nm - 900 nm