



## **Bühler 1027 Stationary easy wall-mount water sampler**

**Product #:** BU1027.XX.XXXXXX  
**AED Price:** Contact Hach

Basic wall-mount sampler using pressure-vacuum technology. With weatherproof plastic housing, suitable for wall-mounting.

The well organised display and numeric keypad enable the sampler to be rapidly programmed. The instrument can be used for time and volume-proportional sampling. Its simple structure means that it requires very little maintenance. It is weatherproof and can be mounted or fixed on walls.

Pressure-vacuum samplers are ISO 5667 compliant and therefore satisfy the requirements for subsequent reproducible analysis.

The samples can be cooled to 4 °C in a stainless steel refrigerator cabinet (optional).

Ideal for straightforward applications involving the sampling of widely different types of wastewater.

**Corrosion-resistant plastic housing**

**Time and volume proportional composite sampling**

**Volume accuracy through pressure-vacuum principle**

**Ideal in combination with stainless steel refrigerator cabinet**

---

### **Specifications**

Cooling system:	External refrigerator option available
Housing material:	Reinforced plastic
Languages user interface:	HR, SI, GR, GB,NL, IT, DK, SE, PL, ES, RU, DE, FR, PT, CZ, HU, RO, SK, BG, TR
Sample container:	1 x 10 L, 1 x 25 L
Sampler type:	Stationary Wall Mount
Warranty:	24 months