



**Be Right™**



## **Moisture analyser IR-35, Denver, 35g capacity, 1mg (0.01%) readability**

<b>Product #:</b>	<b>SMIR35</b>
AED Price:	Contact Hach

with metal rod radiator

Time Preset: 0.1-99 min, Data Interface: RS 232C, Characteristics: Resultdisplay in % Moisture, % Dry Matter, % ATRO, g Remainder, Determinationmode: Automatical Sample Heating: Infraredheating, via metal rod radiator, Results can be printed with opt. printer as a short journal or GLP-protocol in 6 different languages.