



Kit of canisters for 9245-9240 sodium analyser

Product #: 09240=A=8020
AED Price: Contact Hach

Includes:

2 x bottles 0.5L all equipped with caps, stickers and tubing for Reactivtion and AutoCal

1 x Flask PolyEthylene 1L with anti-drop spout

Specifications

Instrument:	9245 & 9240
Measuring range:	N/A
Parameter:	Sodium
Quantity:	N/A
Storage conditions:	Keep containers tightly closed in a dry, cool and well-ventilated place. Protect from moisture.

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)