



Spare parts kit for 2 years, 8810 ISE fluoride analyser

Product #: 368810,93002

AED Price: Contact Hach

Specifications

Composition:

» 1 x 359110,01650 » Stirrer complete
2 x 359090,80000 » Tubing holder cassette
1 x 359110,01400 » Siphon tube for reactor vessel
1 x 359110,01500 » Overflow tube for reactor vessel
2 x 429=110=040 » FKM/FPM O-ring 11x4mm (around above)
2 x 359012,00070 » PTFE washer 18x12.5x3mm (around above)
3 x 359015,20035 » Threaded sleeve PG13,5 Red
1 x 08810=A=0010 » Conn. 90°, G $\frac{1}{4}$ ", Stem 6 mm tube, PP1 x 587=008=002 » Conn. straight, G $\frac{1}{4}$ ", 6/8mm tube, PP1 x 359103,10070 » Conn. 90°, G1/8", 4/6mm tube, PP1 x 359103,10077 » Conn. 90°, G3/8", 10/12mm tube, PP1 x 587=114=003 » Conn. 90°, G3/8", Stem 12mm tube, PP1 x 359103,10132 » Conn. straight, 1/8" MNPT, 6/8mm tube, PP1 x 359110,02300 » Conn. straight, G3/8", Stem 12mm tube, PP10 x 359500,11106 » 'Nipple', 1/16-1/16 thread, stem 1.6mm tube
4 x 359090,70015 » Tube set TYGON grey-grey (12 pieces)
3 x 151399,90002 » Tygon tubing 1,6x3,2mm (per meter)
3 x 359110,04200 » Reagent injection tube, 8810 vessel
3 x 359110,04110 » FEP tubing 1,65x3mm (per meter)
1 x 368483,00000 » Reference electrode 8483-B
1 x 368470,00000 » Fluoride (F-) ISE electrode