



Tensette plus variable electronic pipette

Product #: BBP087
AED Price: Contact Hach

Excellent accurate: Electronic pipette Tensette plus

Electronic Pipette, variable volume 0.2 - 5.0 mL, complete with rack, battery and power supply

Tips are sold separately.

Intuitive menu navigation

Pre-programmed volumes for LCK cuvette tests

Precise and reliable

Easy and safe to handle

Specifications

What's included?: Electronic Pipette 0.2 - 5.0 mL, rack, battery, power supply, basic manual.
Tips are not included and have to be purchased separately.

What's included?

Electronic Pipette 0.2 - 5.0 mL, rack, battery, power supply, basic manual. Tips are not included and have to be purchased separately.