



DRB 200-1 Dry thermostat with 1 block, 9x16mm / 2x20mm

Product #: LTV082.99.30001

AED Price: Contact Hach

Thermostat DRB 200 – the simple solution for all digestions

Pre-programmed for standard digestions (e.g. for COD, P_{tot} or heavy metals) and freely programmable for user specific digestions.

Equipped with a large readable display and a digital timer with an automatic switch-off and acoustic signal. Operational safety guaranteed by protection lid covering the heating block, by the insulated external shell and by an integrated anti-overheating feature. Cuvettes and reaction vessels can be digested at large scale of different temperatures and time settings. Temperature stability and reproducibility ensured.

Great flexibility

Excellent reproducibility

Simple to use

Specifications

| | |
|-------------------------------|--|
| Display resolution: | 5 x 7 dot matrix |
| Display size: | 2 x 16 character |
| Display type: | LCD |
| Heating programmes: | Pre-programmed for 40°C, 100°C, 150°C, 165°C and freely selectable 37-165°C, 1-480 min |
| Heating rate: | From 20 - 150 °C in 10 minutes |
| Heating system: | 1 Heating Block |
| Manual languages: | English, French, German, Italian, Spanish, Dutch, Swedish, Danish, Polish |
| Max. operating humidity: | 90 % |
| Number of cuvettes: | 9 x 16 mm diameter |
| Number of cuvettes (2): | 2 x 20 mm diameter |
| Operating temperature range: | 10 - 45 °C |
| Power requirements (Voltage): | 115/220 V AC |
| Power supply: | 115 V - 230 V +5%/-15% |
| Safety checks: | CE, GS, cTUVus |
| Storage conditions: | -40 - 60 °C |
| Temperature stability: | ± 1 °C in conformity with EN, ISO, EPA methods |
| User interface: | 3 softkeys |
| User programmes: | 6 free temperature/time |

Warranty: 24 months

Weight: 2 kg