



## **DRB 200 Dry thermostat with 2 blocks, 15x16mm, 6x16/4x20mm**

**Product #:** LTV082.99.42001

**AED Price:** Contact Hach

### **Thermostat DRB 200 – the simple solution for all digestions**

Equipped with a large readable display and a digital timer with an automatic switch-off and acoustic signal. Operational safety guaranteed by protection lid covering the heating block, by the insulated external shell and by an integrated anti-overheating feature. Cuvettes and reaction vessels can be digested at large scale of different temperatures and time settings. Temperature stability and reproducibility ensured.

#### **Great flexibility**

Pre-programmed for all

standard digestions and freely programmable for user-specific digestions

#### **Excellent reproducibility**

Temperature stability

#### **Simple handling**

Large, easily readable display and one-key operation for standard

digestions

---

### **Specifications**

Display resolution:	5 x 7 dot matrix
Display size:	2 x 16 character
Display type:	LCD
Heating programmes:	Pre-programmed for 40°C, 100°C, 150°C, 165°C and freely selectable 37-165°C, 1-480 min
Heating rate:	From 20 - 150 °C in 10 minutes
Heating system:	2 Heating Blocks
Manual languages:	English, French, German, Italian, Spanish, Dutch, Swedish, Danish, Polish
Max. operating humidity:	90 %
Number of cuvettes:	4 x 20 mm diameter
Number of cuvettes (2):	21 x 16 mm diameter
Operating temperature range:	10 - 45 °C
Power requirements (Voltage):	115/220 V AC
Power supply:	115 V - 230 V +5%/-15%
Safety checks:	CE, GS, cTUVus

Storage conditions:	-40 - 60 °C
Temperature stability:	± 1 °C in conformity with EN, ISO, EPA methods
User interface:	3 softkeys
User programmes:	6 free temperature/time
Warranty:	24 months
Weight:	2.8 kg