



## LS120 Shaker for surfactant analysis

**Product #:** LQV148.99.10001  
**AED Price:** Contact Hach

For extracting the organic phase, insert for 8 cuvettes, shutoff automatically.

---

### Specifications

|                         |   |
|-------------------------|---|
| Dimensions (H x W x D): | 116 mm x 176 mm x 240 mm  |
| Manual languages:       | English, French, German, Italian, Spanish, Dutch, Swedish, Danish, Polish |
| Parameter:              | Surfactants   |
| Power supply:           | 47 Hz - 63 Hz   |
| Warranty:               | 24 months   |
| Weight:                 | 2.9 kg  |