



Orbisphere luminescent oxygen analyser, panel mount, 85 - 264 V AC

Product #:

41G1-P600

OBSOLETE ITEM

This item is no longer available.

The Orbisphere G1100 oxygen (O_2) luminescent sensor with the Orbisphere 410 offers a new way of monitoring oxygen in high purity water. This on-line process analyser is designed for use in power generation plants running oxygenated treatment (OT).

First luminescent sensor for process ppb level oxygen monitoring

All the benefits of the LDO technology for the power industry

Less than 5 minutes of annual service

Low cost of ownership

Fully automatic in-situ calibration