



Reagent B for TOC process photometer

Product #:

BCF851

Hazardous

AED Price:

Contact Hach



To be ordered as consumable reagent for proper function of the TOC analyser and changed regularly.

Specifications

Special notes:

Reagent B for TOCTAX (5,2 L)

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)