



## Titration Application Pack Water: pH, Alkalinity & Conductivity

**Product #:** AP0002.AT1102  
**AED Price:** Contact Hach  
**Call for ship date**

Application Pack for the determination of pH, Alkalinity and Conductivity in water. Suited for Titralab AT1000 Series, Automatic Titrator Models AT1102/AT1112/AT1122/AT1222.

---

### Specifications

Application:	Municipal Applications
Compatible instruments:	AT1102
	AT1112
	AT1122
	AT1222
Parameter:	pH/Alkalinity & Conductivity
Type:	Application Pack
Warranty:	24 months
What's included?:	Application Pack contains glass syringe, plastic beakers, magnetic stirring bar, pre-programmed application on USB stick, electrode, conical adapter, bottle stopper.

---

### What's included?

Application Pack contains glass syringe, plastic beakers, magnetic stirring bar, pre-programmed application on USB stick, electrode, conical adapter, bottle stopper.