



## Titration Application Pack Water: Ca & Mg Hardness (ISE)

**Product #:** AP0003.AT1112

**AED Price:** Contact Hach

**Call for ship date**

Application Pack for the determination of Calcium and Magnesium Hardness in water. Suited for Titralab AT1000 Series, Automatic Titrator Models AT1112/AT1122/AT1222.

---

### Specifications

**Application:** Municipal Applications

**Compatible instruments:** AT1112

AT1122

AT1222

**Parameter:** Ca & Mg Hardness (ISE)

**Type:** Application Pack

**Warranty:** 24 months

**What's included?:** Application Pack contains glass syringe, plastic beakers, magnetic stirring bar, pre-programmed application on USB stick, electrode, conical adapter, bottle stopper.

---

### What's included?

Application Pack contains glass syringe, plastic beakers, magnetic stirring bar, pre-programmed application on USB stick, electrode, conical adapter, bottle stopper.