



Titration Application Pack Food & Beverage: pH & Total Acidity (Milk, Tomato Sauce, Fruit Juice)

Product #: AP0008.AT1102
AED Price: Contact Hach
Call for ship date

Application Pack for the determination of pH and Total Acidity in food & beverages. Suited for Titralab AT1000 Series, Automatic Titrator Models AT1102/AT1112/AT1122/AT1222.

Specifications

Application:	Food & Beverage Applications
Compatible instruments:	AT1102
	AT1112
	AT1122
	AT1222
Parameter:	pH, Total Acidity
Type:	Application Pack
Warranty:	24 months
What's included?:	Application Pack contains glass syringe, plastic beakers, magnetic stirring bar, pre-programmed application on USB stick, electrode, conical adapter, bottle stopper.

What's included?

Application Pack contains glass syringe, plastic beakers, magnetic stirring bar, pre-programmed application on USB stick, electrode, conical adapter, bottle stopper.