



Titration Application Pack Food & Beverage: Salt

Product #: AP0010.AT1112

AED Price: Contact Hach

Call for ship date

Application Pack for the determination of Salt in food products. Suited for Titrabest AT1000 Series, Automatic Titrator Models AT1112/AT1122/AT1222.

Specifications

Application:	Food & Beverage Applications
Compatible instruments:	AT1112
	AT1122
	AT1222
Parameter:	Salt in food products
Type:	Application Pack
Warranty:	24 months
What's included?:	Application Pack contains glass syringe, plastic beakers, magnetic stirring bar, pre-programmed application on USB stick, electrode, conical adapter, bottle stopper.

What's included?

Application Pack contains glass syringe, plastic beakers, magnetic stirring bar, pre-programmed application on USB stick, electrode, conical adapter, bottle stopper.