



## Titration Application Pack Petrochemistry: TBN (Total Base Number)

<b>Product #:</b>	<b>AP0016.AT1102</b>
AED Price:	Contact Hach
Call for ship date	

Application Pack for the determination of the Total Base Number in Petrochem. Suited for Titrilab AT1000 Series, Automatic Titrator Models AT1102/AT1112/AT1122/AT1222.

---

### Specifications

Application:	Petrochemical Applications
Compatible instruments:	AT1102
	AT1112
	AT1122
	AT1222
Parameter:	TBN (Total Base Number)
Type:	Application Pack
Warranty:	24 months
What's included?:	Application Pack contains glass syringe, plastic beakers, magnetic stirring bar, pre-programmed application on USB stick, electrode, conical adapter, bottle stopper, 1 L glass bottle.

---

### What's included?

Application Pack contains glass syringe, plastic beakers, magnetic stirring bar, pre-programmed application on USB stick, electrode, conical adapter, bottle stopper, 1 L glass bottle.