



## Titration Application Pack Karl Fischer: Moisture Content

**Product #:** AP0014.KF1121  
**AED Price:** Contact Hach

Application Pack for the determination of Moisture content (Karl Fischer). Suited for Titrallab KF1000 Series, Automatic Titrator Model KF1121.

---

### Specifications

Application:	Moisture Content Application
Compatible instruments:	KF1121
Parameter:	Moisture Content (Karl Fischer)
Type:	Application Pack
Warranty:	24 months
What's included?:	Application Pack contains glass syringe, plastic beakers, magnetic stirring bar, pre-programmed application on USB stick, electrode, conical adapter, bottle stopper.

---

### What's included?

Application Pack contains glass syringe, plastic beakers, magnetic stirring bar, pre-programmed application on USB stick, electrode, conical adapter, bottle stopper.