



Total Chlorine Chemkey Reagents (pack of 300)

Product #:	8499400
AED Price:	Contact Hach

Total Chlorine Chemkeys, range 0.04 - 10 mg/L Cl_2 . For SL1000 Portable Parallel Analyser (PPA). 12 boxes of 25 (300 tests total quantity).

Specifications

Instrument:	SL1000
Measuring range:	0.04 - 10 mg/L Cl_2
Method number:	10260
Parameter:	Total Chlorine
Platform:	Chemkeys
Quantity:	pk/300