



## Free Chlorine Chemkey Reagents (box of 25)

|            |              |
|------------|--------------|
| Product #: | 9429000      |
| AED Price: | Contact Hach |

Free Chlorine Chemkeys, range 0.04 - 4.00 mg/L  $\text{Cl}_2$ . For SL1000 Portable Parallel Analyser (PPA).

---

### Specifications

|                  |                                |
|------------------|--------------------------------|
| Measuring range: | 0.04 - 4.00 mg/L $\text{Cl}_2$ |
| Method number:   | 10260                          |
| Parameter:       | Free Chlorine                  |
| Platform:        | Chemkeys                       |
| Quantity:        | pk/25                          |