



## Titration Application Pack Water: Total Hardness (Colour)

**Product #:** AP0004.AT1112

**AED Price:** Contact Hach

Application Pack for the determination of Calcium and Magnesium Hardness in water. Suited for TitraLab AT1000 Series, Automatic Titrator Models AT1112/AT1122/AT1222.

---

### Specifications

Application:	Municipal Applications
Compatible instruments:	AT1112
	AT1122
	AT1222
Parameter:	Total Hardness (Color)
Type:	Application Pack
Warranty:	24 months
What's included?:	Application Pack contains glass syringe, plastic beakers, magnetic stirring bar, pre-programmed application on USB stick, electrode, conical adapter, bottle stopper.

---

### What's included?

Application Pack contains glass syringe, plastic beakers, magnetic stirring bar, pre-programmed application on USB stick, electrode, conical adapter, bottle stopper.