



Titration Application Pack Petrochemistry: R-SH (Thiol)

Product #: AP0017.AT1112

AED Price: Contact Hach

Application Pack for the determination of R-SH (Thiol) in Petrochem. Suited for Titralab AT1000 Series, Automatic Titrator Models AT1112/AT1122/AT1222.

Specifications

Application:	Petrochemical Applications
Compatible instruments:	AT1112
	AT1122
	AT1222
Parameter:	R-SH (Thiol)
Type:	Application Pack
Warranty:	24 months
What's included?:	Application Pack contains glass syringe, plastic beakers, magnetic stirring bar, pre-programmed application on USB stick, electrode, conical adapter, bottle stopper, 1 L glass bottle.

What's included?

Application Pack contains glass syringe, plastic beakers, magnetic stirring bar, pre-programmed application on USB stick, electrode, conical adapter, bottle stopper, 1 L glass bottle.