



Nitrate Standard Solution, 1 mg/L, 500 mL

Product #: 204649
AED Price: Contact Hach
Available

Ion Selective Electrode (ISE) standard solutions for laboratory analyses

ISE Standard Solution. Nitrate, 1.00 mg/L ± 0.01 mg/L $\text{NO}_3\text{-N}$ (NIST). 500 mL bottle

Traceable to NIST standard reference materials

Always performing measurement results you can trust is essential. The quality of the standard solutions that are used are key.

Design for use with Hach meters portfolio and other meters.

The basis for accurate ISE measurements is good care and calibration of the Ion Selective Electrode.

Reliable packaging

Supplied in an airtight packaging ensuring guaranteed shelf-life.

Specifications

| | |
|---------------------|--|
| Concentration: | 4.43 mg/L ± 0.04 mg/L NO_3 |
| Description: | Nitrogen-Nitrate Standard Solution as $\text{NO}_3\text{-N}$ (NIST) |
| Package Type: | Bottle |
| Parameter: | Nitrate |
| Platform: | Solution |
| Quantity: | 500 mL |
| Special remarks: | Observed shelf life after opening, following recommended storage and use: 3 months |
| Storage conditions: | 10 - 25 °C (protect from light) |