



Quality control standard metals for drinking water, low range, 500 mL

Product #:	2833749
AED Price:	Contact Hach
Available	

Ready to use - no dilutions required. Mixed standard includes: Copper: 1 mg/L Cu, Iron: 0.3 mg/L Fe, Manganese: 0.1 mg/L mg/lMn. 500 mL Bottle.

Specifications

Concentration:	Iron: 0.3 mg/L Fe
Description:	Metals Drinking Water (Low Range) Mixed Parameter QC Standard
Platform:	Solution