



## Filter, glass fiber, pore size 2.7 µm, diameter 70 mm, 100/pk

|                   |                |
|-------------------|----------------|
| <b>Product #:</b> | <b>2164300</b> |
| AED Price:        | Contact Hach   |
| Available         |                |

Made from 100% borosilicate glass fibers. Best used in general laboratory work and as a prefilter for membrane filtration.

---

### Specifications

|            |                  |
|------------|------------------|
| Diameter:  | 70 mm            |
| Material:  | Glass Microfiber |
| Pore size: | 2.7 µm           |
| Quantity:  | 100 /pk          |