



## QL filters, pore size 8 - 12 $\mu\text{m}$ , diameter 90 mm, circular, 100/ pk, qualitative

<b>Product #:</b>	<b>50655</b>
AED Price:	Contact Hach
Available	

QL Filters are recommended for use in various environmental and wet chemistry reference methods. Extremely pure cellulose filters of nearly 100% alpha content are specially processed to yield negligible extractable levels in the determination and identification of unknown sample materials.

---

### Specifications

Diameter:	90 mm
Material:	Cellulose
Quantity:	100 /pk