



## Flask, volumetric, polypropylene, 250 mL

Product #: 1406046

**OBSOLETE ITEM**  
This item is no longer available.

Polypropylene flask and screw cap. Durable with excellent chemical resistance, with the same accuracy as glass. Permanent color graduation line with marking spot. Calibrated “to contain/deliver” at 20°C. The 200/203 capacity flask has dual marks useful in DO determination via Winkler Titration.

---

### Specifications

Capacity:	250 mL
Certifications:	Not Applicable
Housing material:	Polypropylene
Units:	1 each