



Beaker, Griffin, Low Form, Glass, 600 mL

Product #:	50052
AED Price:	Contact Hach
Available	

Griffin style (low form) beakers are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing. Bull nose lip and molded pour spout minimise spillage in liquid transfer. Autoclavable, non-toxic and biologically inert. Volumes and graduations are approximate.

Specifications

Capacity:	600 mL
Housing material:	Glass
Material:	Glass, Griffin, Low Form
Units:	each